REF 11751

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1.22 SpO2 Testing Summary

SpO2 accuracy, and low perfusion testing were conducted using the Masimo SET Technology and sensors with Masimo optical sources and receivers. The Iradimed Corporation fiberoptic MRI SpO2 sensors utilize the same Masimo active optical components. Masimo Rainbow SET technology (using sensors with Masimo optical sources and receivers) has been validated for no motion accuracy in human blood studies on healthy adult volunteers in induced hypoxia studies in the range of 70-100% SpO2 against a laboratory CO-Oximeter and ECG monitor. This variation equals plus or minus one standard deviation which encompasses 68% of the population.

Masimo Rainbow SET technology (using sensors with Masimo optical sources and receivers) has been validated for motion accuracy in human blood studies on healthy adult volunteers in induced hypoxia studies while performing rubbing and tapping motions, at 2 to 4 Hz at an amplitude of 1 to 2 cm and a nonrepetitive motion between 1 to 5 Hz at an amplitude of 2 to 3 cm in induced hypoxia studies in the range of 70-100% SpO2 against a laboratory co-oximeter and ECG monitor. This variation equals plus or minus one standard deviation which encompasses 68% of the population.

Masimo Rainbow SET technology (using sensors with Masimo optical sources and receivers) has been validated for low perfusion accuracy in bench top testing against a Biotek Index 2 simulator and Masimo's simulator with signal strengths of greater than 0.02% and a % transmission of greater than 5% for saturations ranging from 70 to 100%. This variation equals plus or minus one standard deviation which encompasses 68% of the population.

1.23 Specifications.

Oxygen Saturation Display Range: 0 to 100% SpO₂

Pulse Rate Display Range: 30 to 240 beats per minute (BPM)

SpO₂ Accuracy (Arms)a:

FO Sensor 70-100% ± 2 digits

Pulse Rate Accuracy (Arms)a:

No Motion ±3 digits, 30-240 BPM

Low Perfusion ±3 digits, 30-240 BPM

a.) ±1 Arms represents approximately 68% of Measurements

SpO₂/Pulse Rate Response Time: < 10 seconds

Measurement Wavelengths and Output Power:b

Red: 660 nm @ 0.8 mW maximum average

Infrared: 910 nm @ 1.2 mW maximum average

b.) This information is especially useful for clinicians performing photodynamic therapy.

SpO₂ Temperature (Operating): 0o to +40o C (32o F to 104o F)

Temperature (Storage/Transportation): -30o to +50o C (-22o F to 122o F)

SpO₂ Humidity (Operating): 10 to 90% noncondensing

Humidity (Storage/Transportation): 10 to 95% noncondensing

SpO₂ Altitude (Operating): Up to 12,000 meters (40,000 feet)

SpO₂ Sensor Accuracy

SpO₂: 70 to 100% ±2 digits (Arms*).

Pulse Rate Accuracy: ±3 digits (BPM)

*±1 Arms represents approximately 68% of measurements.

Probe Compliance:

This product complies with ISO 10993-1.

1.24 The Masimo Set® Principles of Operation.

The Masimo SET® MS pulse oximeter is based on three principles:

1. Oxyhemoglobin and deoxyhemoglobin differ in their absorption of red and infrared light (spectrophotometry).
2. The volume of arterial blood in tissue and the light absorbed by the blood changes during the pulse (plethysmography).
3. Arterio-venous shunting is highly variable and that fluctuating absorbance by venous blood is a major component of noise during the pulse.

The Masimo SET® MS pulse oximeter as well as traditional pulse oximetry determines SpO₂ by passing red and infrared light into a capillary bed and measuring changes in light absorption during the pulsatile cycle. Red and infrared light-emitting diodes (LEDs) in oximetry sensors serve as the light sources, a photodiode serves as the photodetector.

Traditional pulse oximetry assumes that all pulsations in the light absorbance signal are caused by oscillations in the arterial blood volume. This assumes that the blood flow in the region of the sensor passes entirely through the capillary bed rather than through any arterio-venous shunts. The traditional pulse oximeter calculates the ratio of pulsatile absorbance (AC) to the mean absorbance (DC) at each of two wavelengths, 660 nm and 905 nm:

$$S(660) = AC(660)/DC(660)$$

$$S(905) = AC(905)/DC(905)$$

The oximeter then calculates the ratio of these two arterial pulse-added absorbance signals:

$$R = S(660)/S(905)$$

This value of R is used to find the saturation SpO₂ in a look-up table built into the oximeter's software. The values in the look-up table are based upon human blood studies against a laboratory co-oximeter on healthy adult volunteers in induced hypoxia studies.

The Masimo SET® MS board pulse oximeter assumes that arterio-venous shunting is highly variable and that fluctuating absorbance by venous blood is the major component of noise during the pulse. MS board decomposes S(660) and S(905) into an arterial signal plus a noise component and calculates the ratio of the arterial signals without the noise:

$$S(660) = S1 + N1$$

$$S(905) = S2 + N2$$

$$R = S1/S2$$

Again, R is the ratio of two arterial pulse-added absorbance signals and its value is used to find the saturation SpO₂ in an empirically derived equation into the oximeter's software. The values in the empirically derived equation are based upon human blood studies against a laboratory co-oximeter on healthy adult volunteers in induced hypoxia studies.

The above equations are combined and a noise reference (N') is determined:

$$N' = S(660) - S(905) \times R$$

If there is no noise N' = 0: then S(660) = S(905) x R which is the same relationship for the traditional pulse oximeter.

The equation for the noise reference is based on the value of R, the value being sought to determine the SpO₂. The MS board software sweeps through possible values of R that correspond to SpO₂ values between 1% and 100% and generates an N' value for each of these R-values. The S(660) and S(905) signals are processed with each possible N' noise reference through an adaptive correlation canceler (ACC) which yields an output power for each possible value of R (i.e., each possible SpO₂ from 1% to 100%). The result is a Discrete Saturation Transform (DST™) plot of relative output power versus possible SpO₂ value as shown in the following figure where R corresponds to SpO₂ = 97%:


Figure 1-10. Saturation Charting.

The DST plot has two peaks: the peak corresponding to the higher saturation is selected as the SpO₂ value. This entire sequence is repeated once every two seconds on the most recent four seconds of raw data. The MS board SpO₂ therefore corresponds to a running average of arterial hemoglobin saturation that is updated every two seconds.

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APPENDIX A SPECIFICATIONS

2460/2465 GENERAL SYSTEM REQUIREMENTS	
General Features	Pulse Oximeter for use in the MRI
ELECTRICAL CHARACTERISTICS	
HI/LO Line Voltage Requirements	100 to 240 VAC +/-10%, 50/60 Hz
Power Sources Available	Internal battery power with separate AC charger/power supply
Power Consumption	< 5 Volt-Amperes @ 120 VAC nominal (<72 VA maximum during charging)
Battery Type	Rechargeable Lithium Polymer Pack, 14.8 v at 5.8 Ah
Battery Capacity (2460 Only)	> 20 Hours
Battery Charge Time	< 9 hours to 95% capacity
Battery Cycle Life	> 200 charge/discharge cycles
Patient Leakage Current (2460 Only)	< 20 uA RMS
Chassis Leakage Current	< 300 uA RMS
Chassis Ground Impedance to Earth Ground	< 0.1 ohms (with power supply)
MECHANICAL CHARACTERISTICS	
Dimensions D x W x H	3 x 6 x 9 Inches (7.6 x 15.2 x 22.9 cm)
Total Weight	2.7 LBS (1.2 Kg); 5 LBS (2.3 Kg) with Battery.
Operating Temperature Range	+5° to +40° C
Storage Temperature Range	-40° TO +70° C
Operating Relative Humidity Range	0% to 80% RH, non-condensing
Storage Relative Humidity Range	0% to 95% RH, non-condensing
Pole Clamp Mounting (Pole Diameter) Range (2460 Only)	1 inch to 1.5 inch (25 to 38 mm) diameter

PERFORMANCE CHARACTERISTICS	
MRI Performance	
MRI Magnet Compatibility (2460 Monitor)	0.2 to 3.0 Tesla MRI Systems
MRI Magnet Compatibility (2465 Remote)	None Specified - 2465 Remote Unit is NOT intended for use inside the Magnet Room.
MRI Larmor Frequency Interference	No interference at compatible Larmor Frequencies (8.4 through 128 MHz)
Magnetic Field Limit	30,000 Gauss (3 Tesla magnetic field line) - Minimal ferrous material used inside Oximeter (< 15 grams).
PRODUCT STANDARDS REQUIREMENTS	
STANDARDS COMPLIANCE: IEC 60601-1	Yes
STANDARDS COMPLIANCE: IEC 60601-1-2	Yes
STANDARDS COMPLIANCE: UL 2601	Yes
STANDARDS COMPLIANCE: ISO 9919	Yes
REMOTE COMMUNICATIONS	
FCC Certification	Part 15.247, no license required
ETS (European) Certification	brETS 300.328, no license required
Rated RF Output Power	+18 dBm (default), 24 dBm (maximum) Service Selectable.
Frequency Range	2431.3 to 2474.5 MHz (Iradimed Default) Avoids 802.11 bands 1 and 2 France: 2400 to 2450 MHz with Radio programming selection. Service Selectable.
Number of Channels	85 U.S. Channels. Note: A limited frequency mode is available for France. Service Selectable.
OTHER	
Replacement Fuse Type (1120 Only)	1 Ampere, 3AG, 250V Slo-Blo or equivalent.

APPENDIX B REPAIR

All repairs on products under warranty must be performed by Iradimed Corporation personnel, or an authorized Iradimed Corporation Service and Repair Center. Unauthorized repairs will void the warranty.

If a Oximeter fails to function properly or requires maintenance, contact Iradimed Corporation Technical Service at 1-407-677-8022 within the U.S., +001-407-677-8022 from outside the U.S. (during normal business hours EST), or by E-mail at techsupport@iradimed.com. Iradimed Corporation Technical Service will advise you of the corrective action required. If you are advised to return the Oximeter to Iradimed Corporation for repair, please do the following:

- a. Obtain a Return Authorization Number. This will ensure proper routing and facilitate timely repair of your Oximeter.
- b. Clean Oximeter prior to shipment. Do not ship contaminated product to IRadimed Corporation for repair.
- c. Package the Oximeter with adequate protection. If available, use the original carton and packing materials in which the Oximeter was shipped from Iradimed Corporation.
- d. Include a brief description of the problem as well as the name, address and phone number of the person to be contacted for additional information.
- e. Include a purchase order with the Oximeter being returned if it is out of warranty; Iradimed Corporation Technical Services can advise you of your Oximeter's warranty status if need be. Repairs will be made at Iradimed Corporation's current list price for the replacement part(s) plus a reasonable labor charge.
- f. Ship the Oximeter transportation prepaid, to the location specified by your Iradimed Corporation Service Representative with the Return Authorization Number written on the outside of the shipping carton. Repairs will be made, normally, within two (2) weeks and the Oximeter will be returned to you prepaid.

To ensure full reliability, it is recommended that all repairs be made by an Iradimed Corporation Authorized Service and Repair center. For repair at your facility, a competent individual experienced in the repair of Oximeters can repair the Oximeter only IF it is authorized by Iradimed Corporation Technical Service prior to the repair.

CAUTION: No repair should ever be attempted by anyone not having a complete knowledge of the repair of Iradimed Corporation Oximeters. Only replace damaged parts with components manufactured or sold by Iradimed Corporation. Contact the Iradimed Corporation Technical Service and Repair Center for service and technical assistance.

APPENDIX C

WARRANTY INFORMATION

Iradimed Corporation warrants this product, other than its expendable parts, to be free from defects in materials and workmanship for a period of twelve (12) months from the date of original delivery to the buyer or to buyer's order, provided that same is properly operated under conditions of normal use, and that periodic maintenance and service is performed. This same warranty is made for a period of thirty (30) days on expendable parts. This warranty will become null and void if product has been repaired other than by Iradimed Corporation, or if the product has been subject to misuse, accident, negligence or abuse.

Iradimed Corporation's sole obligation under this warranty is limited to repairing a product which has been reported to Iradimed Corporation's Technical Service Center during normal business hours and shipped transportation prepaid. Iradimed Corporation is not liable for any damages including, but not limited to, incidental damages, consequential damages or special damages.

This warranty is in lieu of any other warranties, guarantees or conditions, including merchantability or fitness for a particular purpose. The remedies under this warranty are exclusive and Iradimed Corporation neither assumes nor authorizes anyone to assume for it any other obligation in connection with the sale or repair of its products.

IRADIMED CORPORATION PRODUCTS CONTAIN PROPRIETARY COPY WRITTEN MATERIAL; ALL RIGHTS ARE RESERVED BY IRADIMED CORPORATION.

APPENDIX D
MANUFACTURERS TECHNICAL DECLARATION

EMC Information Tables as required by EN 60601-1-2:2000 Clause 6

In accordance with EN 60601-1-2:2000 Medical Electrical Equipment – Part 1-2: General requirements for safety – Collateral standard: Electromagnetic Compatibility – Requirements and tests

1. “Medical Electrical Equipment needs special precautions regarding EMC and needs to be installed and put into service according to the EMC information provided in the Accompanying Documents” (the following tables).
2. “Portable and Mobile RF Communications Equipment can affect Medical Electrical Equipment”.
3. “The equipment or system should not be used adjacent to or stacked with other equipment and that if adjacent or stacked use is necessary, the equipment or system should be observed to verify normal operation in the configuration in which it is used”.

The following tables (as required in EN 60601-1-2:2000) provide information regarding the Electromagnetic Compatibility (EMC) of this product and its accessories.

Table 201—Guidance and manufacturer’s declaration— electromagnetic emissions—for all EQUIPMENT and SYSTEMS

Guidance and manufacturer’s declaration—electromagnetic emissions		
The Imagox 2460 System is intended for use in the electromagnetic environment specified below. The customer or the user of the Imagox 2460 System should assure that it is used in such an environment.		
Emissions Test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	The Imagox 2460 System must emit electromagnetic energy in order to perform its intended function (remote communications within a specific band for WLAN; i.e. 2.431 to 2.474 GHz). Nearby electronic equipment may be affected.
RF emissions CISPR 11	Class B	The Imagox 2460 System is suitable for use in all establishments other than domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Not applicable	
Voltage fluctuations/flicker emissions IEC 61000-3-3	Not applicable	

Imagox 2460 System

Table 202—Guidance and manufacturer's declaration— electromagnetic immunity—for all EQUIPMENT and SYSTEMS

Guidance and manufacturer's declaration—electromagnetic immunity			
The Imagox 2460 System is intended for use in the electromagnetic environment customer or the user of the Imagox 2460 System should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±6 kV contact ±8 kV air	±6 kV contact ±8 kV air	Floors should be wood, concrete, or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/burst IEC 61000-4-4	±2 kV for power supply lines ±1 kV for input/output lines	±2 kV for power supply lines ±1 kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±1 kV differential mode ±2 kV common mode	±1 kV differential mode ±2 kV common mode	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions, and voltage variations on power supply input lines IEC 61000-4-11	<5 % UT (>95 % dip in UT) for 0.5 cycle 40 % UT (60 % dip in UT) for 5 cycles 70 % UT (30 % dip in UT) for 25 cycles <5 % UT(>95 % dip in UT) for 5 sec	<5 % UT (>95 % dip in UT) for 0.5 cycle 40 % UT (60 % dip in UT) for 5 cycles 70 % UT (30 % dip in UT) for 25 cycles <5 % UT(>95 % dip in UT) for 5 sec	Mains power quality should be that of a typical commercial or hospital environment. If the user of the Imagox 2460 System requires continued operation during power mains interruptions, it is recommended that the Imagox 2460 System be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

NOTE—UT is the a.c. mains voltage prior to application of the test level.

Table 203—Guidance and manufacturer’s declaration— electromagnetic immunity—for LIFE-SUPPORTING EQUIPMENT and SYSTEMS

Guidance and manufacturer’s declaration—electromagnetic immunity			
The Imagox 2460 System is intended for use in the electromagnetic environment customer or the user of the Imagox 2460 System should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz outside ISM bands ^a	3 V	Portable and mobile RF communications equipment should be used no closer to any part of the Imagox 2460 System, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance $d = 1.17 \sqrt{P}$
Radiated RF IEC 61000-4-3	10 Vrms 150 kHz to 80 MHz outside ISM bands ^a	10 V	$d = 1.20 \sqrt{P}$
	10 Vrms 80 MHz to 2.5 GHz	10 V/m	$d = 1.20 \sqrt{P}$ 80 MHz to 800 MHz $d = 2.33 \sqrt{P}$ 800 MHz to 2.5 GHz where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). ^b Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, ^c should be less than the compliance level in each frequency range. ^d Interference may occur in the vicinity of equipment marked with the following symbol: IEC 60417, No. 417-IEC5140, "Source of Non-Ionizing Radiation" Symbol
NOTE 1—At 80 MHz and 800 MHz, the higher frequency range applies.			
NOTE 2—These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.			
<p>^a The ISM (industrial, scientific, and medical) bands between 150 kHz and 80 MHz are 6.765 MHz to 6.795 MHz; 13.553 MHz to 13.567 MHz; 26.957 MHz to 27.283 MHz; and 40.66 MHz to 40.70 MHz</p> <p>^b The compliance levels in the ISM frequency bands between 150 kHz and 80 MHz and in the frequency range 80 MHz to 2.5 GHz are intended to decrease the likelihood that mobile/portable communications equipment could cause interference if it is inadvertently brought into patient areas. For this reason, an additional factor of 10/3 is used in calculating the recommended separation distance for transmitters in these frequency ranges.</p> <p>^c Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast, and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the MRidium 3850 System is used exceeds the applicable RF compliance level above, the MRidium 3850 System should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the MRidium 3850 System.</p> <p>^d Over the frequency range 150 kHz to 80 MHz, field strengths should be less than [V/m]</p>			

Imagox 2460 System

Table 205—Recommended separation distances between portable and mobile RF communications equipment and the EQUIPMENT or SYSTEM—for LIFE-SUPPORTING EQUIPMENT and SYSTEMS

Recommended separation distances between portable and mobile RF communications equipment and the Imagox 2460 System				
The Imagox 2460 System is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the Imagox 2460 System can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the Imagox 2460 System as recommended below, according to the maximum output power of the communications equipment.				
Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m			
	150 kHz to 80 MHz outside ISM bands $d = 1.17 \sqrt{P}$	150 kHz to 80 MHz in ISM bands $d = 1.2 \sqrt{P}$	80 MHz to 800 MHz $d = 1.2 \sqrt{P}$	800 MHz to 2.5 GHz $d = 2.33 \sqrt{P}$
0.01	0.12	0.12	0.12	0.23
0.1	0.38	0.38	0.38	0.77
1	1.17	1.2	1.17	2.33
10	3.8	3.8	3.8	7.67
100	11.67	12	11.67	23.33

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be determined using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1—At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2—The ISM (industrial, scientific, and medical) bands between 150 kHz and 80 MHz are 6.765 MHz to 6.795 MHz; 13.553 MHz to 13.567 MHz; 26.957 MHz to 27.283 MHz; and 40.66 MHz to 40.70 MHz.

NOTE 3—An additional factor of 10/3 is used in calculating the recommended separation distance for transmitters in the ISM frequency bands between 150 kHz and 80 MHz and in the frequency range 80 MHz to 2.5 GHz to decrease the likelihood that mobile/portable communications equipment could cause interference if it is inadvertently brought into patient areas.

NOTE 4—These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.

DECLARATION OF CONFORMITY

Application of Council Directive(s) R&TTE Directive 1999/5/EC

Regulation(s) to which Conformity is Declared R&TTE Directive 1999/5/EC

Standard(s) to which Conformity is Declared EN 300 328, IEC 60601-1, IEC 60601-1-2

Manufacturer's Name Iradimed Corporation

Manufacturer's Address 7457 Aloma Ave., Winter Park, FL 32792 USA

Importer's Name Refer to accompanying Packing Slip

Importer's Address Refer to accompanying Packing Slip

Type of Equipment Imagox MRI Pulse Oximeter System and accessories

Model No(s). 2460 Series and 2465 Series

Serial No(s). Refer to accompanying Packing Slip

Year of Manufacture Refer to accompanying Packing Slip

Certification Method(s) Self-certified

I, the undersigned, hereby declare that the equipment specified above conforms to the above Directive(s), Regulation(s) and/or Standard(s).

Place Winter Park, FL

Date October 31, 2011

(Signature)

Francis Casey

(Full Name)

Regulatory Affairs

(Position)

**Declaration of Conformity with Regard to the EU Directive 1999/5/EC
(R&TTE Directive)**

Scope

The information in this document is applicable to the Imagox 2460 MRI Oximeter System using the 2460 MRI Pulse Oximeter with 2.4 GHz FHSS Transceiver, and 2465 Remote Display Controller with the 2.4 GHz FHSS Transceiver, and other related products as they come to market.

The equipment operates in the bands (or portions of) 2400 to 2483.5 MHz. National regulations may require that operations be limited to a portion of the frequency ranges identified above. See the National Restrictions section for full details.

Declaration of Conformity

Czech (CS) Česky	Toto zařízení je v souladu se základními požadavky a ostatními odpovídajícími ustanoveními směrnice 1999/5/ES
Danish (DA) Dansk	Dette udstyr er i overensstemmelse med de væsentlige krav og andre relevante bestemmelser i Direktiv 1999/5/EF.
Dutch (NL) Nederlands	Dit apparaat voldoet aan de essentiële eisen en andere van toepassing zijnde bepalingen van de Richtlijn 1999/5/EG
Estonian (ET) Eesti	See seade vastab direktiivi 1999/5/EÜ olulistele nõuetele ja teistele asjakohastele säteletele.
English (EN)	This equipment is in compliance with the essential requirements and other relevant provisions of Directive 1999/5/EC.
German (DE) Deutsch	Dieses Gerät entspricht den grundlegenden Anforderungen und den weiteren entsprechenden Vorgaben der Richtlinie 1999/5/EG
Greek (EL) Ελληνική	Αυτός ο εξοπλισμός είναι σε συμμόρφωση με τις ουσιαστικές απαιτήσεις και άλλες σχετικές διατάξεις της Οδηγίας 1999/5/ΕΚ.
Finnish (FI) Suomi	Tämä laite täyttää direktiivin 1999/5/EY olennaiset vaatimukset ja on siinä asetettujen muiden laitetta koskevien määräysten mukainen.
French (FR) Français	Cet appareil est conforme aux exigences essentielles et aux autres dispositions pertinentes de la Directive 1999/5/CE
Hungarian (HU) Magyar	Ez a készülék teljesíti az alapvető követelményeket és más 1999/5/EK irányelvben meghatározott vonatkozó rendelkezéseket.
Icelandic (IS) Íslenska	Þetta tæki er samkvæmt grunnkröfum og öðrum viðeigandi ákvæðum Tíðskipunar 1999/5/EB
Italian (IT) Italiano	Questo apparato è conforme ai requisiti essenziali ed agli altri principi sanciti dalla Direttiva 1999/5/CE
Latvian (LV) Latviski	Šī iekārta atbilst Direktīvas 1999/5/EK būtiskajām prasībām un citiem ar to saistītajiem noteikumiem.
Lithuanian (LT) Lietuvių	Šis įrenginys tenkina 1999/5/EB Direktyvos esminius reikalavimus ir kitas šios direktyvos nuostatas.

Maltese (MT) Malti	Dan l-apparat huwa konformi mal-htigiet essenzjali u l-provedimenti l-ohra rilevanti tad- Direttiva 1999/5/KE
Norwegian (NO) Norsk	Dette utstyret er i samsvar med de grunnleggende krav og andre relevante bestemmelser i EU-direktiv.1999/5/EC
Polish (PL) Polski	Urządzenie jest zgodne z ogólnymi wymaganiami oraz szczególnymi warunkami określonymi Dyrektywą UE: 1999/5/WE
Portuguese (PT) Português	Este equipamento está em conformidade com os requisitos essenciais e outras provisos relevantes da Directiva 1999/5/CE

Applicable Standards

The following standards were applied during the assessment of the product against the requirements of the Directive 1999/5/EC:

- Radio: EN 300 328
- EMC: EN 60601-1-2
- Safety: EN 60601-1-1 and EN 60601-2-24

CE Marking

For the 2.4 GHz wireless products such as the Imagox 2460 System using the 2460 MRI Pulse Oximeter with 2.4 GHz FHSS Transceiver, and 2465 Remote Display Controller with the 2.4 GHz FHSS Transceiver, the following CE mark and class-2 identifier are added to the equipment:

National Restrictions

In the majority of the EU and other European countries, the 2.4 GHz band has been made available for the use of ISM (Industrial, Scientific, and Medical) communication. Table 1 provides an overview of the regulatory requirements in general that are applicable for the 2.4 GHz band.

Later in this section you will find an overview of countries in which additional restrictions or requirements or both are applicable.

The requirements for any country may evolve. Iradimed Corporation recommends that you check with local authorities for the latest status of their national regulations for 2.4 GHz ISM (Industrial, Scientific, and Medical) communication.

Frequency Band: 2400 - 2483.5 MHz

Max Power Level (EIRP): 100 mW

Indoor and Outdoor Use: Yes

The following sections identify countries having additional requirements or restrictions than those listed above:

France

For 2.4 GHz, the output power is restricted to 10 mW EIRP when the product is used outdoors in the band 2454 - 2483.5 MHz. There are no restrictions when used in other parts of the 2.4 GHz band. For use in France, channel selection is limited to 2400 to 2450 MHz. Check <http://www.arcep.fr/> for more details.

Italy

This product meets the National Radio Interference specifications and the requirements in the National Frequency Allocation Table for Italy. Operating wireless equipment requires a "general authorization" unless it is operated within the boundaries of the owner's property. Please check with <http://comunicazioni.it/it> for more information.

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Latvia

Any outdoor usage of the 2.4 GHz band requires an authorization form the Electronic Communications Office. Please check <http://www.esd.lv> for more information.

Antennas

The 2.4 GHz wireless ISM (Industrial, Scientific, and Medical) communication products described in this document have dedicated antennas which should not be removed or replaced with any antenna that is not supplied directly from Iradimed Corporation or its representatives.

Operating Frequency

The operating frequency is pre-set at the factory during manufacturing. Manual setting of the operating frequency can only be done by qualified service personnel with a service

password. Instructions are provided in the Service documentation regarding the appropriate choices and settings for specific geographic areas. It is important that the operating frequency is correctly configured to meet the local regulations.

Changing Output Power

The output power is pre-set to 63 mW at the factory during manufacturing. Manual setting of the output power above 100 mW can only be done by qualified service personnel with a service password. The only available output power choices are Level 1 (10 mW), Level 2 (63 mW – factory default), or Level 3 (250 mW). Instructions are provided in the Service documentation regarding the appropriate choices and settings for specific geographic areas. It is important that the output power is correctly configured to meet the local regulations.

Note If you still have questions regarding the compliance of these products, or you cannot find the information you are looking for on the Iradimed Corporation website, <http://iradimed.com>, please send an email request to Iradimed Corporation at:

techsupport@iradimed.com

APPENDIX E ACCESSORIES

Accessories:

- 1119 - MRI IV Pole (Non-Magnetic) I.V. Pole, 3.2 CM Dia. (1.25 in.)
- 1120 - MRI Power Supply
- 1121 - AC Power Adapter Power Cord, 1 Meter (3 ft.) (length), medical grade.
- 1122 - AC Adapter/To Oximeter Interconnect Cable, 3 Meters (10 ft.)(length).
- 1128 - AC Power Cord, 3.05 m (length), Medical Grade
- 1133R - MRI-Compatible Replacement Battery Pack
- 1175 - Imagox™ 2460 Series Operators Manual
- 1176 - Imagox™ 2460 Series Service Manual CD
- 1170 - Fiberoptic SpO₂ Sensor Kit
- 1171 - Finger Grip Kit Includes:
 - 1171A - Adult Figer Grip
 - 1171P - Infant/Pediatric Grip
 - 1171L - Adult Large Finger Grip
- 1172 - Fiberoptic SpO₂ Cable Pole Mount
- AM10 - Imagox™ 2460 Pulse Oximeter System Software Update Card

WARNING: Use only listed accessories with the Imagox™ Model 2460 MRI Pulse Oximeter System.

APPENDIX F

1119 I.V. POLE ASSEMBLIES AND PARTS DESCRIPTIONS

General Precautions

- Federal law in the USA restricts this device to sale by or on the order of a physician.
- For safe operation, use only Iradimed Corporation recommended MRI-compatible or MRI-safe accessories.
- Always secure the IV Pole wheel locks after positioning within the MRI Magnet Room.
- Product damage may occur unless proper care is exercised during unpacking and installation. Use only recommended assemblies and parts. Any replacement component must be non-magnetic for safe operation.
- Refer all service to Iradimed Corporation, or Authorized Service Representatives.
- A maximum number of two (2) Oximeters can be safely used on any one (1) IV pole, mounted no more than 54 inches [137cm] from the floor.

NOTE: There are several 1119 I.V. Pole configurations. When using the following assembly instructions, please match the appropriate instructions to your specific 1119 I.V. Pole.

Routine Maintenance

Periodically check all mounting hardware. Tighten as necessary for optimal operation.

Cleaning the Mounting Assembly

1. The mounting assembly may be cleaned with most mild, non-abrasive solutions commonly used in the hospital environment (e.g. diluted bleach, ammonia, or alcohol solutions).
2. The surface finish will be permanently damaged by strong chemicals and solvents such as acetone and trichloroethylene.
3. Do not use steel wool or other abrasive material to clean the mounting assembly.
4. Damage caused by the use of unapproved substances or processes will not be covered by warranty. We recommended that you test any cleaning solution on a small area of the mounting assembly that is not visible to verify compatibility.
5. Never submerge the roll stand or allow liquids to enter the mounting assemblies. Wipe any cleaning agents off the mounting assembly immediately, using a water-dampened cloth. Dry all mounting assemblies thoroughly after cleaning.

CAUTION: Iradimed Corporation makes no claims regarding the efficacy of the listed chemicals or processes as a means for controlling infection. Consult your hospital's infection control officer or epidemiologist. To clean or sterilize mounted instruments or accessory equipment, refer to the specific instructions delivered with those products.


IV POLE CONFIGURATION
1119 I.V. POLE ASSEMBLY AND PARTS DESCRIPTIONS

NOTE: The following tools are required for the assembly of this IV Pole:


Small Screws 1/8" HEX

Big Screws 3/16" HEX

Bolthead 1/2" Wrench

STEP1: ASSEMBLE THE BOTTOM POST+ BUSHING ASSEMBLY ONTO THE ROLLSTAND BASE


STEP2: SCREW THE HOOKS INTO THE SWIVEL PIECE, SECURE THE SWIVEL TO THE TOP POST USING THE 3x FLAT HEAD SCREWS PROVIDED. NOTE THAT THE TRI-HOLE PATTERN IS ASYMETRICAL, AND THEREBY ONLY ONE CORRECT HOLE-ALIGNMENT CONFIGURATION IS POSSIBLE.


PRIVATE AND CONFIDENTIAL
 THE REPAIR INSTRUCTIONS OF AMCO ACCESSORIES AND EQUIPMENT VEHICLE REPAIRS ARE SUBJECT TO SEVERAL PATENT RIGHTS AND TRADE MARKS. THE INFORMATION CONTAINED HEREIN IS FOR YOUR INFORMATION ONLY AND IS NOT TO BE REPRODUCED OR TRANSMITTED IN ANY FORM OR BY ANY MEANS, ELECTRONIC OR MECHANICAL, INCLUDING PHOTOCOPYING, RECORDING, OR BY ANY INFORMATION STORAGE AND RETRIEVAL SYSTEM, WITHOUT PERMISSION IN WRITING FROM AMCO ACCESSORIES.
 ALL RIGHTS TO OTHER TRADEMARKS ARE RESERVED.

9 9

		SHEET NO. 00000 PRINTED IN CANADA 11/03/19/10
TITLE: ROLLSTAND CUSTOMER MANUAL		
PART: STEP 1 & 2		
DATE: 3/1	REV: 001/19/10	REV: 001/19/10
ISSUED TO: 3/1	ISSUED TO: 3/1	ISSUED TO: 3/1

STEP3: SLIDE THE TOP WEIGHT ON THE POST AND REST IT ON BASE ASSEMBLY


STEP4: ATTACH THE BOTTOM WEIGHT TO THE TOP WEIGHT USING SX FLAT-HEAD SCREWS


PRIVATE AND CONFIDENTIAL
 THIS PRINT IS PROPERTY OF AUSTON ATTORNEYS
 AND IS LOANED IN CONFIDENTIAL SUBJECT TO RETURN
 UNREPRODUCED AND WITH THE AUSTON ATTORNEYS THAT
 NO TOYONS KEYS TO BE MADE WITHOUT THE CONSENT
 OF AUSTON ATTORNEYS
 ALL RIGHTS TO BE OBSERVED UNDER ART. 17(2)(b)

REFERENCE: ISSUANCE CUSTOMER MANUAL	
PART TO: STEP 3 & 4	
DRAWN BY: EC	DATE: 04/01/2010
CHECKED BY:	DATE:
00019-001	00019-001


NOTE:
TIGHTEN THE BOLT UNTIL THE
STICKER DISSAPPEARS
UNDER THE WEIGHT

STEP5: SECURE THE POST ASSEMBLY TO THE
WEIGHT ASSEMBLY USING THE LARGE AND SMALL
WASHER AND THE BOLT, IN THE SEQUENCE
SHOWN IN THE ASSEMBLY

STEP5: SLIDE THE TOP POST ASSEMBLY
ONTO THE CONNECTOR OF THE
BOTTOM ASSEMBLY AND LOCK THE
PIECES TOGETHER USING THE
HOLE+PLUNGER MECHANISM
NOTE: PUSH THE PLUNGER PIN WITH
A TOOL

PLUNGER PIN
(ENSURE TIP OF PLUNGER
IS PROTRUDING FROM
THE POST)

PRIVATE AND CONFIDENTIAL
THE FRONT STRUT ROD IS A RIVET ACCESSORY
AND A STEERING COMPONENT. IT IS NOT TO BE
LINED UP WITH THE FRONT STRUT ROD AND
NO COPIES ARE TO BE MADE WITHOUT THE CONSENT
OF AMCO ACCESSORIES

SYSTEM OF CARE SAFETY SYSTEM SAFETY SYSTEM SAFETY SYSTEM	
DOCUMENT RADIUM CUSTOMER MANUAL	STEP 5 & 6
PART NO. DATE: MAR/19/21C	PART NO. DATE: MAR/19/21C
ALL RIGHTS RESERVED OR VARIATIONS ARE RESERVED	100319-PC1